

## Technical Specifications for Harmony LITE

Electrical	100-120 VAC, 20A, 50/60 Hz 220-240 VAC, 10A, 50/60 Hz
Weight	63 Kg (139 lbs)
Dimensions	35 cm wide x 40 cm deep x 110 cm high (14" x 16" x 43.3")

Module	Applications	Technology	Wavelength (nm)	Pulse Width	Spot Size (cm <sup>2</sup> )	Pulse Repetition Rate (Hz)	Energy Density (J/cm <sup>2</sup> )
High Power UV	Pigment Restoration	Pulsed UV	300-380	30,40, 50 msec	6.4	1/2	200-5000 mJ/cm <sup>2</sup>
Acne Pro	Acne Clearance	AFT	420-950	30,40, 50 msec	6.4	1/2	5-25
Dye-VL	Vascular & Pigmented lesions	AFT	500-600	10, 12, 15 msec	3	1/2	Up to 15
VL/PL	Vascular & Pigmented lesions	AFT	540-950	10, 12, 15 msec	3	2/3	5-30
SSR Pro	Skin Rejuvenation	AFT	570-950	10, 12, 15 msec	3	2/3	5-30
SHR Pro	Super hair removal	In-Motion AFT	NIR	8 msec	5	3	1-7
NIR	Muscles & joints pain relief	Pulsed NIR	NIR	NA	6.4	NA	5-25
SST Pro	Scar revision & Treatment of striae	In-Motion AFT	NIR	4 msec	6.4	5	1-3.5
LED	Various skin conditions	LED	Yellow	Pulsed / CW	NA	NA	NA

\* The parameters are subject to change

### About Alma Lasers

Alma Lasers, Ltd. is a global developer, manufacturer and provider of laser, light-based, radiofrequency and ultrasound devices for aesthetic and medical applications. Since 1980, the founders of Alma Lasers have been at the forefront of innovative multi-technology / multi-application systems designed to meet the unique needs of today's practitioners.

Alma's mission is to provide modular, cost-effective and high-performance systems that enable practitioners to confidently offer safe, effective and profitable aesthetic treatments to their patients.

**Innovative Technology**  
Safe, effective treatments with the industry's most dynamic product portfolio, distinguished by expandable modular platforms that adapt to your changing needs and future breakthroughs.

**Innovative Marketing**  
Consumer-outreach focus, marketing assistance kits, A.L.M.A. marketing consulting support and Office by Alma online practice management store & resource center, all coordinated to attract new patients and increase ROI.

**Innovative Learning**  
An unparalleled commitment to customer care and knowledge, demonstrated through in-office training, hands-on workshops, intensive webinars, advanced clinical training and our online education resource center.

**Innovative Support**  
On-call service technicians standing by 24/7 and world-class Depot Service with guaranteed overnight delivery of loaner equipment to eliminate downtime. \*certain exclusions apply

© 2012 Alma Lasers Ltd. All Rights Reserved. Alma Lasers, its logos, Harmony, and Pain-Free Hair-Free are the trademarks or registered trademarks of Alma Lasers Ltd. in the United States and/or other countries. Product specifications are subject to change without notice.

**Alma Lasers®**  
Wellbeing Through Technology®

**US Headquarters**  
485 Half Day Road  
Suite 100  
Buffalo Grove, IL 60089  
Tel +1-224-377-2000  
Fax +1-224-377-2050  
marketing@almalasers.com

**International Headquarters**  
14 Halamish St.  
Caesarea Industrial Park  
Caesarea, Israel 38900  
Tel +972-4-627-5357  
Fax +972-4-627-5368  
info@almalasers.com

PBHA10011204



# Better Investment, Better Technology

The Harmony<sup>LITE</sup> is the world's newest expandable IPL platform for aesthetic light treatments. Its modular design provides your practice multiple distinct cosmetic solutions all built into one compact unit. With the Harmony<sup>LITE</sup>, you can purchase only the IPL handpieces you require today and add new technologies as your practice grows.

The Harmony<sup>LITE</sup> gives you ultimate flexibility to treat a wide variety of patient-requested procedures. The system's modular handpieces are designed to deliver effective and reliable results with minimal to no patient discomfort and downtime.

Harmony<sup>LITE</sup> provides the largest selection of IPL applications for safe and effective treatments.

- ▶ Extensive range of IPL wavelengths, many with integrated contact cooling
- ▶ High output energy for greater efficacy
- ▶ High patient and clinical staff satisfaction
- ▶ Easily transportable
- ▶ Customised to fit your practice
- ▶ Affordable with a high return on investment

▶ **IN-Motion™**

- High Power UV
- Acne Pro
- Dye-VL
- VL/PL
- SSR Pro
- SHR Pro
- SST Pro
- NIR
- LED



## Results You Can Rely On



Before

After

**Acne Pro**

Photos courtesy: Elizabeth VanderVeer, MD



Before

After

**Dye-VL**

Photos courtesy: Clinical department of Alma Lasers Ltd.



Before

After

**SSR Pro**

Photos courtesy: Fernando Stengel, MD



Before

After

**SHR Pro**

Photos courtesy: Mira Volvovsky, MD